

# Nanopore Conductance Analysis GUI — User Guide

Last updated: 2025-08

---

## 1) What this app does

A desktop GUI for analyzing I-V data from nanochannels/nanopores. It lets you: - Load I-V curves (multiple files). - Convert units consistently (original → target) and re-scale plots/data in place. - Compute conductance (global and user-selected voltage ranges). - Detect the zero-current crossing ( $E_{\text{total}}$  : left/mid/right) and annotate plots. - Switch between **IV** and **G-V** views with one click. - Compare experimental vs theoretical conductance/ conductivity ( $\sigma$ -C). - Calculate drift-diffusion quantities, redox/Em, and ion mobilities ( $\mu^+/\mu^-$ ,  $\mu^+$ ,  $\mu^-$ ) for monovalent and general salts. - Compute enclosed area of a looped I-V curve and log frequency-normalized metrics. - Manage an **Electrolyte details** table (e.g., KCl/LiCl conductivity vs concentration) for quick reference in calculations. - Export figures and data to PNG/CSV.

**Non-blocking plotting:** Live plotting (if enabled) runs decoupled from acquisition in a separate module/process. Interacting with the plot window (drag/resize/zoom) will **not** affect measurement timing.

**Communication latency note:** Due to VISA/SCPI round-trip and instrument measurement time, even with *dwell per step* = 0, the effective time spacing between points is typically ~25 ms.

---

## 2) Requirements

- **Python** 3.9+
  - **Libraries:** `tkinter`, `numpy`, `pandas`, `matplotlib`, `mplcursors`
  - Optional (not required for analysis): `pyvisa` (for instrument discovery in other tooling)
- 

## 3) Launching the app

Run the main analysis script (e.g., `python analysis_gui.py`). The window has: - **Left panel:** file controls, unit conversion, analysis actions, export buttons. - **Right panel:** active plot **or** results table (context-dependent).

---

## 4) Supported data formats

- **Excel** (`.xlsx`): first two columns are **Voltage (V)** and **Current (A)**.
- **CSV** (`.csv`): header row is skipped; column 2 → **Voltage**, column 3 → **Current**.

- **Text** ( `.txt` ): two numeric columns → **Voltage, Current** (whitespace-separated).

Files are cleaned for `NaN` rows before use. If a file has fewer than two numeric columns, you'll be asked to fix the input.

## 5) Quick start (typical workflow)

1. **Load files** → multiple I-V curves can be added/removed at any time.
2. **Set units** → choose the **original units** of your data (e.g., V, mV, A, nA) and your **target units**. Click **Rescale** to convert in place. (This respects the chosen *original* unit as the baseline—no hidden resets.)
3. **Inspect IV plot** → the app automatically:
  4. Fits a global linear conductance (shows slope on the figure).
  5. Finds zero-current crossings and labels **Ettotal (left, mid, right)**.
  6. Displays hover readouts via `mplcursors`.
7. **Switch mode** → use the top-right buttons **IV / G-V**.
8. **IV mode**: current vs voltage, colored by sweep direction; legends auto-update.
9. **G-V mode**: conductance envelopes vs voltage; y-axis is **Conductance (S)**.
10. **Selected-range conductance** → click **Calculate Conductance in Selected Region**.
11. A small window opens (separate from the main UI). Enter your voltage range.
12. The app computes the local linear conductance and shows a zoomed plot.
13. Results are written into an editable table and can be exported.
14. **Enclosed area** (optional) → **Calculate Enclosed Area** opens a plot and an editable table for polygon area of looped I-V curves. You can append the corresponding frequency for normalization.
15. **Export** → save current figure as PNG and data as CSV. "Save all" writes every active series/table.

## 6) Zero-current crossings (Ettotal) and table behavior

- The midpoint crossing is auto-filled in the table after loading a curve.
- Selecting **Ettotal (left/mid/right)** from the dropdown updates the **currently selected row** in the drift-diffusion table (it does **not** add a new row).
- These values can be edited later and are used directly in mobility and Em workflows.

## 7) Drift-diffusion & mobility module

This module computes Em, mobilities, and ratios based on your inputs.

**Inputs** - Concentrations on each side ( $C_H$ ,  $C_L$ ) and, if needed, activity coefficients  $\gamma$ . - Temperature  $T$ . - Valences  $z^+$ ,  $z^-$  (monovalent or general salts). - Either **Em directly** or a **redox potential** (from which Em is derived).

**Outputs** -  $\mu^+/\mu^-$ ,  $\mu^+$ ,  $\mu^-$  and conductivity  $\sigma$ , using the Henderson form (monovalent) or the general-salt formulation.

**Usage tips** - For KCl/LiCl, you can enter **activity coefficients (y)** per side if required. - All results are placed into an editable table. Use **Update** to recompute for a row after editing inputs.

---

## 8) Electrolyte details (conductivity table)

Click **Enter Electrolyte Details** to open a side table with prefilled KCl/LiCl conductivity (S/m) at common concentrations. You can: - Edit values and add new salts/concentrations. - Save the table for later runs. - Use these values directly in theoretical conductance/ $\sigma$  computations.

---

## 9) Reference window

Open the **References** window for a concise explanation of the underlying formulas and clickable links to source papers/handbooks. This is a separate window you can open anytime; reading it will not interrupt ongoing analysis.

---

## 10) Unit conversion

Use **Attach Unit Conversion UI** to pick both **original units** and **target units**. Press **Rescale** to convert all loaded data and any dependent tables/plots. This applies to IV plots, G-V, and  $\sigma$ -C figures.

---

## 11) Exports & data management

- **Save Figure** → PNG of the current plot.
  - **Export Data** → CSV for the visible table or the active curves.
  - **Save All** → batch export of all displayed curves/tables.
  - Tables are editable: double-click to change a cell; use **Save** to persist.
- 

## 12) Performance & timing

- **Non-blocking UI:** Live plots run in a separate module/process; UI actions won't block acquisition.
  - **25 ms floor:** Even with dwell per step = 0, communication + measurement overhead limits the effective point interval to  $\approx 25$  ms.
-

## 13) Troubleshooting

**"File must contain at least two numeric columns."** Check the file format and column order (V first, I second). For CSV, ensure the first line is a header (it will be skipped).

**"Plot shows wrong units."** Confirm you set **original** units correctly before clicking **Rescale**.

**"Selected-range conductance not saved."** Ensure the sub-window says "Result written to table." If you closed it before saving, re-run the calculation.

**"G-V looks empty."** If you switched from IV after hiding data, click **IV** then **G-V** again to refresh envelopes.

**Performance feels slow.** Close unused figures; avoid extremely large CSV/Excel files; confirm you aren't plotting millions of points at once.

---

## 14) Tips

- Use hover (mouse over points) to see precise coordinates.
  - Keep a consistent naming scheme for files so exported CSVs map back to raw data easily.
  - For batch comparisons, load files together, rescale once, then export "All."
- 

## 15) Known limits

- CSV dialects with irregular delimiters may fail auto-sniffing; open and save as standard CSV if needed.
  - Extremely noisy I-V data may degrade linear fits. Use **Selected-range conductance** to focus on regions with good linearity.
- 

## 16) Change log (highlights)

- Live plot moved to its own module/process (non-blocking).
  - Etotal selection updates the current table row (no extra rows created).
  - Added **Selected-range conductance** window and zoomed preview.
  - Added **Electrolyte details** table and  $\sigma$ -C integration.
  - Unified unit-conversion UI across IV, G-V, and  $\sigma$ -C.
- 

**Questions or requests?** Feel free to annotate this guide with your notes; we'll keep it updated as the app evolves.